

A.R.C. LASER

HYBRID SYSTEMS FOR ENT

The all new A.R.C. hybrid laser system

consisting of a green 532 nm Laser „NUVOLAS“ and the CO₂ body “C-LAS” as an option you easily can add a Scanner system „WEASEL” in addition.



GREEN

NUVOLAS green 532 nm laser

CO₂ LASER

with articulated arm and various hand pieces



optional:

SCANNER

high speed ultra fine scanner „Weasel“



optional:

High precision micro-manipulator with a working distance up to 500 mm



Perform your surgery with the micromanipulator in „Line of Sight”

